Description du dispositif

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Tel: +90 (232) 469 14 44

www.wingmed.com.tr - info@wingmed.com.tr

Dogam Medikal company manufacturer of othopedic equipment under Wingmed brand declares that, marked with the CE sign is Class I medical device, of rule 1, consistent with the requirements of Regulation (EU) 2017/745 of the European Parliament and of the Council of 5 April 2017 on medical devices. Conformity assessment was carried out in accordance with the requirements of Regulation (EU) 2017/745 of the European Parliament and of the Council.

tr	Medikal Cihaz	Tek hasta Çoklu Kullanım	Talimat kitapçığını okuyun
en	Medical Device	Single patient Multiple Use	Read the instruction leaflet
fr	Dispositif médical	Un seul patient à usage multiple	Lire attentivement la notice
ru	Медицинский прибор	Многократное использование для одного пациента	Внимательно прочтите инструкцию

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Device Description

This device is a type of orthosis designed to protect and support the ligaments in the knee joint. The device is designed only for current indications and accurate body measurement.

Device Content

Textile Components: Polyester-Polyamide-Elastane

Rigid Components: Silicone-Steel

Device Application

Insert your leg through the device, centering the patella and positioning the inner label towards the top of the knee.

Indications

- Increased intra-articular fluid
- Side ligament damage and kneecap pain
- In rehabilitation after surgical intervention

Contraindication

- Do not use the product in cases where the diagnosis is not certain.
- Do not contact the product directly with the injured skin.
- Do not use in case of allergy to any of the components of the product.
- Do not use in venous and lymphatic conditions.

Important Information & Warnings

Check the integrity of the device before each use. Do not use if the device is damaged. It is recommended that a healthcare professional supervise. Do not use the Device while driving. This device is only suitable for use by one patient and should not be used by another person. In case of discomfort, pain, abnormal volume and discoloration of the limb, remove the device and consult a healthcare professional. The duration of daily use and application should be indicated by the health care professional. Substances containing oil should not be brought into contact with (ointments, creams, lotions).

Device Maintenance & Storage

It should be washed in warm water and hands with mild detergent. To preserve the long-term adhesive strength of the velcro tapes, the velcro tapes should be closed before washing. It should be dried in the open air without hanging and should not be exposed to heat (sun, heating) when drying. Do not dry clean. Keep the device clean in its original box when not in use.

S	31-34
M	34-37
L	37-40
XL	40-44
XXL	44-48

Size

DİZ ÇEVRESİ
KNEE CIRCUMFERENCE
CIRCONFÉRENCE DU GENOU
ОКРУЖНОСТЬ КОЛЕНА (cm)